

RNA
Catalog #RN-00-001



RNA

Instructions for Use

Catalog #RN-01-001

Use Restrictions

This product is distributed for research use only. The use of this product is not approved for human or veterinary use. Do not use this product with in vitro diagnostics.

Buyer has no rights to transfer products, components, or materials made using these products, or to use these products for commercial purposes. Commercial purposes include: use of products or their components in manufacturing; use of products or their data components to provide a service, information or data; use of products or their components for therapeutic or diagnostic purposes; resale of products or their components.

Safety

This product may contain human source material. Treat as potentially infectious. Handle at Biological Safety Level 2 to minimize exposure to potentially infectious agents.

User Information

The RNA is isolated from Cook MyoSite's Human Skeletal Muscle Derived Cells (SkMDC), **Cook MyoSite Catalog #SK-1111**. The RNA provides the researcher with access to RNA of SkMDC populations to save valuable time and cost. The RNA is suspended in RNase-free water at 0.5 µg/µL. The RNA has been DNase-treated to ensure minimal DNA contamination. This product has been quantified and determined to be of high quality and ready to use in any downstream application such as real-time RT-PCR and array analysis.

Additional Information

Primary cultures have a finite lifespan and limited number of population doublings in vitro. Cook MyoSite cryopreserved cells and associated reagents are tested prior to shipment for contaminating agents. Contamination of cell cultures may affect cell growth, function and behavior.

For detailed information concerning Quality Control testing and specifications, please refer to Certificate of Analysis.

COOK MYOSITE INC.
105 Delta Drive
Pittsburgh, PA 15238
412.963.7380